

AquaSafe

WSL25 Plus

Water Safety

Laboratory

Instruction Manual

Contents

Section 1 : Introduction

Section 2 : Kit Contents

Section 3 : WSL25 Incubator Instructions

Section 4 : Membrane Filtration Instructions

Section 5 : Turbidity Tube Instructions

Section 6 : Comparator Instructions

Section 7 : pH Pocket Tester Instructions

Section 8 : EC Pocket Tester Instructions

Introduction

The Aquasafe WSL25 Plus is a fully portable field laboratory for the detection of microbiological and physico-chemical water quality parameters.

The lab is housed in a waterproof wheeled carry case with telescopic handle and comprises the following :

- Integrated single chamber digital incubator complete with all accessories and consumables for carrying out 250 tests for Faecal and Total Coliforms.
- Comparator complete with accessories and consumables for 250 tests of Free & Total Chlorine.
- Pocket meter for pH
- Pocket meter for EC
- Turbidity Tube for turbidity testing in the range 5 to 500NTU.



Kit Contents



Fig 1.0 Main Case components





Fig 1.1 Tray components

Main Case Components

Item	Qty	Description
1	1	Turbidity Tube 5-500NTU
2	1	Digital Incubator
3	1	Sampling Cable with Carabiner
4	1	Pen
5	1	Multi-Purpose Screwdriver
6	1	Digital Thermometer
7	1	Vacuum Pump
8	5	Lauryl Sulphate Broth Sachet
9	1	Silicone grease
10	1	Calibration Pack
11	1	Vehicle Charging Cable
12	5	Sample Collection Bag (with dechlorination tablet
13	1	Box Membrane Filters (100 Pack)
14	1	Phenol Red reagents (50 Pack)
15	1	DPD No.1 reagents (50 Pack)
16	1	DPD No.3 reagents (50 Pack)
17	2	Comparator Tube
18	1	HydroLite HL101 pH Pocket Tester
19	1	HydroLite HL102 EC Pocket Tester
20	1	Removable Tray

Tray Components

Item	Qty	Description
1	5	Sterilised water for Membrane Lauryl Sulphate Broth preparation
2	1	5ml Sterile syringe
3	1	30ml Methanol
4	1	Membrane Filtration Unit including, waste beaker, sample beaker, measuring funnel, sintered glass disk & 3 silicone gaskets
5	3	Container of Pads (100 Pack)
6	1	Forceps
7	1	Vacuum Tube
8	1	Eye Glass
9	1	Pad Dispenser
10	1	Comparator



Section 3

AquaSafe

MSI25

Microbiological

Safety

Incubator

Instruction Manual

1.0 Introduction

The Aquasafe MSI25 is a portable incubator for the incubation of microbiological samples prepared using the membrane filtration method. The incubator is primarily designed to be used with the supplied 54mm x 3.5mm aluminium petri dishes which are suitable for 47mm membrane filter pads, but the incubator can also be used with 55mm pre-prepared plastic petri-dishes.

The incubator is supplied with 20 Petri dishes as standard but has a capacity for up to 25 dishes.

The incubator has the option to run 37°C, 44°C and user defined temperature profiles for periods of 1 – 24 hours.



2.0 Incubator set-up

The incubator is integrated into the waterproof carry case with a rechargeable battery and power supply.

The power pack is comprised of a 12Volt 12Ah Sealed **Non-Spillable** Lead acid battery. The battery must be electrically isolated during transport on-board aircraft so to facilitate this the power supply is disconnected from the battery by removal of the positive(red) terminal connection and the DC Jack is removed from the incubator socket and fitted with a cap. Prior to use the connector labelled A must be connected to the positive terminal of the battery labelled B (accessed by removing the white panel) and the DC Jack (labelled as such) should be plugged into the socket on the side of the incubator after removal of the cap.

Please note when transporting via aircraft these connections should be removed to prevent any risk of short circuit during transit.

It is recommended to charge the battery fully prior to use. To do this remove the power supply unit from the kit without overstressing the cable and plug into a 240V AC outlet. The indicator on the charger will illuminate red, orange or green. The states indicate battery condition as follows :

Red : discharged.

Orange : Partially charged.

Green : Fully charged.

To charge the battery from empty takes approximately 4-6 hours.

If desired the incubator can be used whilst the battery is being charged however the battery will charge at a lower rate.

The incubator is also supplied with a vehicle cigarette lighter cable. To use this the incubator can be removed from the kit by lifting it up and unplugging the DC jack on the right hand side of the unit. The cigarette lighter cable is then plugged into this socket. 12V or 24V vehicle systems may be used.

3.0 Operating instructions

3.1 Incubator controls

Incubator controls

The incubator has an on/off switch and power connector on the right hand side. The on/off switch is pressed in to turn on and depressed to turn off.

The incubator can be run from a power supply from 12 – 24Volts.

The user interface is via the three switches and the LCD display on the top of the incubator.



Petri-Lok Cassette.

The Petri-Lok cassette consists of a rack and a spring loaded top to hold the petri-dishes in place. To remove the cassette, twist the cassette to the left and lift out. The top of the cassette is removed by pushing down on the top whilst holding the rack, twisting to the left and then removing the top. To re-fit align the arrow with either of the two slots in the top of the rack, push down and twist lock into place.

3.2 Getting started

It is assumed that the user is familiar with the preparation of samples using the membrane filtration method and as such this manual does not cover this however, instructions on membrane filtration are available on request.

The user should familiarise themselves fully with the set-up, control and calibration of the incubator prior to carrying out tests for the first time.

Turn on the power to the incubator using the switch located on the right hand side just above the 12V DC power connector. This is accessed by the recess in the foam to the right of the incubator.

The incubator will beep and the screen will turn on.

The incubator may be in one of three operating modes on power up :

- 1 : New run
- 2: Partial incubation cycle (idle)
- 3: Partial incubation cycle (Run)

If the incubator is in mode 1 or 2 the screen shown in Fig1.0 will be displayed. If the incubator is in mode 3 the incubator will return to the last know state and auto start the incubation cycle from where it was before the power was removed. Mode 3 is primarily intended for unattended power failure conditions.

If Fig 1.0 is displayed and the incubator was in mode 2 then the previous incubation cycle can be restarted from where it was left by pressing restart. This will be discussed later.


The screen will be displayed for 60seconds before defaulting to the restart condition.

Pressing the right hand switch will reset the incubator ready to set-up and start a new incubation cycle and the screen will change to the one shown in Fig 1.1.


Pressing the left hand button will step through the screens as illustrated in the cyclic sequence shown in Figures 1.1 to 1.8.

Each screen and options will now be discussed in further detail.

3.3 Menu sequence



3.3.1 Start screen (Figure 1.1)



This screen shows the current status of the incubator.

In the state shown the incubator is in idle mode and not running.

3.3.2 37°C Program (Fig 1.2)

20.9°C	24:00
->37°C	MEM 10.5v

This mode is the pre-programmed default 37°C incubation cycle. In this mode the incubator will run an 18 hour incubation cycle at 37°C. The 18 hour timer can be altered to any setting between 0 and 24 hours. Adjusting the incubation time is described in section 3.3.6

Pressing and holding the middle switch for 5 seconds will start the incubation cycle and the screen will change to the following :

*20.9°C	24:00
->37°C	RUN 10.5v

The asterisk preceding the temperature will switch on and off indicating the status of the heater.

To stop or pause the incubation cycle press and hold the middle switch again for 5 seconds. The display will return to the one shown in Fig1.2.

3.3.3 44°C Program (Fig1.3)

20.9°C		24:00
->44°C	MEM	10.5v

This mode is the pre-programmed default 44°C incubation cycle. In this mode the incubator will run an 18 hour incubation cycle at 44°C. The 18 hour timer can be altered to any setting between 0 and 24 hours. Adjusting the incubation time is described in section 3.3.6

Pressing and holding the middle switch for 5 seconds will start the incubation cycle and the screen will change to the following :

*20.9°C		24:00
->44°C	RUN	10.5v

The asterisk preceding the temperature will switch on and off indicating the status of the heater.

To stop or pause the incubation cycle press and hold the middle switch again for 5 seconds. The display will return to the one shown in Fig1.2.

3.3.4 User defined Program (Fig 1.4)

20.9°C		24:00
U->##°C	MEM	10.5v

This mode will run an incubation cycle at the user defined temperature for 18 hours. Adjusting the user defined temperature is described in section 3.3.7. The 18 hour timer can be altered to any setting between 0 and 24 hours. Adjusting the incubation time is described in section 3.3.6

Pressing and holding the middle switch for 5 seconds will start the incubation cycle and the screen will change to the following :

*20.9°C		24:00
U->##°	RUN	10.5v

The asterisk preceding the temperature will switch on and off indicating the status of the heater.

To stop or pause the incubation cycle press and hold the middle switch again for 5 seconds. The display will return to the one shown in Fig1.2.

3.3.5 Current Run status (Fig 1.5)

20.9°C	24:00
New	<SEL 10.5v

In this mode pressing the middle switch will display the following screen which will flash intermittently :

15:45	TIMER A
RST	SAVE EXIT

The time remaining is displayed in the top left corner. This can be reset to the programmed default cycle time by pressing the left switch. After the switch press it will jump to the screen shown in Fig1.6. Pressing save will save the current time and jump to the screen shown in Fig1.6. Exiting will return from this mode.

3.3.6 Incubation cycle time setting (Fig 1.6)

20.9°C	24:00
Tim	<ADJ 10.5v

In this mode the default cycle time can be changed from 18 hours to anywhere between 0 and 24 hours.

To adjust the time press the middle switch. The screen will change to :

18:00	TIMER A
-	Save +

Pressing the left switch will decrease the time and pressing the right switch will increase the time. When the desired time has been set press and hold the middle switch to save the settings.

3.3.7 User defined temperature setting (Fig1.7)

20.9°C	24:00
U##	<ADJ 10.5v

In this mode the user defined incubation temperature can be adjusted to a temperature between 20°C and 50°C.

To adjust the temperature press the middle switch. The screen will change to :

##.°C	User Temp
-	Save +

Use the left switch to decrease the temperature and the right switch to increase the temperature. When the desired temperature is set press and hold the middle switch to save the settings.

3.3.8 Calibration (Fig 1.8)

20.9°C		24:00
Cal	<ADJ	10.5v

In this mode the incubator is calibrated.

To calibrate the unit the following equipment is required:

- Digital thermometer with 100mm x 4mm stainless steel probe (Not included)
- Trace2o Incubator calibration pack (Included)
- Petri-Lok petri-dish cassette.
- 10 Aluminium Petri dishes.

Place 10 empty Petri dishes into the bottom of the cassette.

Then place the calibration pack into the cassette.

Fit the spring loaded top and insert the thermometer through the hole into the top until it is fully inserted. See Fig 3.1.

Fit the cassette into the incubator chamber.

Set the incubator to the desired mode of operation – 37/44/User, and start the incubation cycle. Wait at least 1 hour.



Fig 3.1

Compare the displayed temperature on the incubator screen to that shown on the thermometer. If there is a difference > 0.2°C the calibration should be adjusted as follows :

Exit the run mode and step through the screens until the calibration screen is displayed. Press the middle switch to enter calibration mode. The following screen will be displayed.

##.°C	Calibrate
-	Save +

Adjust the temperature up or down using the left and right switches until it matches the temperature displayed on the thermometer. Press the middle switch to save the settings.

Step through the screens until the mode previously been used is displayed and start the incubation cycle again.

Allow the incubator time to adjust the temperature which could take up to 30mins and check the calibration again. If the original temperature was significantly different to the thermometer then the calibration steps may need to be repeated two or three times.

4.0 Care and Maintenance

The Aquasafe incubator is designed to require minimal maintenance although from time to time cleaning will be required.

To clean the incubator outer surfaces and plastics use a cloth moistened with a mild soap solution or an alcohol impregnated cloth or wipe. **Do not** use acetone or other solvents.

The incubator is splash proof but it is not fully waterproof so it should not be left outside in wet conditions and **must not** be submersed in water.

The incubator chamber is sealed at the bottom so spills will not affect the unit but should be cleaned as soon as possible. The Aluminium chamber can be cleaned with Alcohol, soap or solvent based cleaners.

The Petri-Lok rack can be steam sterilised if required but it is not recommended to steam sterilise the top.

5.0 Guarantee and Assistance

Trace2o hope that the Aquasafe incubator will give many years of trouble free operation but in the event of a technical problem occurring the Aquasafe incubator is covered by Trace2o Ltds standard Guarantee terms and conditions available via email or via download from www.trace2o.com.

In the event that any technical assistance is required Trace2o Customer service department will be happy to assist. Contact details as follows :

Trace2o Ltd
The Technology Centre
Wagtech Court
Station Road
Thatcham
Berkshire. RG19 4HZ
UK.
T: 01635 866772
E: Technical@Trace2o.com

Section 4

AquaSafe

Membrane Filtration

Instruction Manual

Contents

Contents	2
SECTION 1: ASSEMBLY OF FILTRATION APPARATUS	3
SECTION 2: PREPARING BACTERIOLOGICAL MEDIA IN A LABORATORY FACILITY	5
SECTION 3: PREPARING BACTERIOLOGICAL MEDIA IN THE FIELD	6
SECTION 4: SAMPLING	7
SECTION 5: USE OF BACTERIOLOGICAL MEDIA	9
SECTION 6: ASEPTIC PROCEDURES	10
SECTION 7: PROCESSING SAMPLES FOR COLIFORM ANALYSIS	12
SECTION 8: COUNTING COLIFORMS AND RECORDING THE RESULT	16
SECTION 9: SELECTING THE OPTIMUM VOLUMES FOR MEMBRANE FILTRATION	17

SECTION 1: ASSEMBLY OF FILTRATION APPARATUS



1. Hand Vacuum Pump
2. Filtrate Flask / Waste Beaker
3. Sampling Cup and cable
4. Graduated Aluminium Funnel
5. Membrane Support & Holder
6. Sealing Gaskets
7. Glass Sintered Disc (Membrane Support)

Place one gasket within the recess of the support holder and press into place. Place the sintered disc, smooth side upmost, into the centre of the gasket. Push the other two gaskets into place around the sintered support disc.



The gaskets, disc and membrane holder

The graduated funnel screws clockwise into position. Do not over-tighten the funnel.



Fitting the graduated funnel

SECTION 2: PREPARING BACTERIOLOGICAL MEDIA IN A LABORATORY FACILITY

Membrane Lauryl Sulfate Broth: For 50 tests, dissolve 7.62g of Membrane Lauryl Sulfate Broth (one sachet) in 100mL deionised water.

The broth is supplied in a pre-weighed sterile sachet, with indicating silica gel that will turn from orange to green, to indicate moisture penetration.

Tear open the sachet and remove the silica gel with the forceps. If the silica gel is green, discard the sachet without using.

Pour the entire contents of the sachet into 100mL of deionised water.

Gently heat the mixture to ensure that the powder is fully dissolved, but do not boil.

After ensuring that the 125mL plastic bottles provided contain no residues of previous MLSB or cleansing agent, pour the prepared medium carefully into the bottles.

Replace bottle lids but leave them slightly loose - do not fully tighten.

Sterilise in an autoclave at 121°C for 10 minutes, or place bottles in a pressure cooker and maintain in steam at pressure for 15 minutes. Remove the bottles, allow them to cool to room temperature, fully tighten the tops and then store in a cool, dark place.

When the media has cooled to room temperature, pour sufficient MLSB (2mL) onto each membrane pad to saturate the pad.

When the pad is fully saturated, decant any excess MLSB as waste

Media Ampoules: Media Ampoules are pre-sterilised ampoules containing 2 mL of dissolved media. They have the advantage of convenience and of always being sterile. These ampoules are available for Faecal Coliform Counts (pack of 50) and Total Coliform Counts (pack of 50). Simply unscrew the cap, pour the media onto the pad and discard the empty ampoule.

SECTION 3: PREPARING BACTERIOLOGICAL MEDIA IN THE FIELD

Membrane Lauryl Sulfate Broth: For 50 tests, dissolve 7.62g of Membrane Lauryl Sulfate Broth (one sachet) in 100mL deionised water.

The broth is supplied in a pre-weighed sterile sachet, with indicating silica gel that will turn from orange to green, to indicate moisture penetration.

Tear open the sachet and remove the silica gel with the forceps. If the silica gel is green, discard the sachet without using.

Open a bottle of 'Sterilised water for broth preparation' and carefully pour the entire contents of the sachet into the bottle.

Replace the bottle lid tightly, and vigorously shake the mixture to ensure that the powder is fully dissolved.

Use the syringe provided to dispense 2mL of prepared broth onto each membrane pad.

SECTION 4: SAMPLING

Rivers and streams

Take the sample as near as possible to the fastest flow – this will typically be found towards the centre of the body of water. Avoid taking samples from too close to the bank, where the water may be still and unrepresentative.

Care must be taken not to introduce floating debris, or solid matter from the banks of the water course, into the water sample. Therefore, it may be preferable to attach the sampling cable to the sterilised sampling cup and take the sample from a bridge or other overhanging location. Alternatively, the cup may be cast out into the water from the edge and pulled slowly and carefully back towards the operator.

Chlorinated drinking water

When storing samples of water in bottles for analysis at a later date, from sources that contain residual chlorine, such as treated drinking water, the bottled sample must be dechlorinated, i.e. the residual chlorine has to be removed to prevent further chlorination (and killing of bacteria) whilst in transit.

To dechlorinate the water, the supplied sample collection bags contain dechlorination tablets.

N.B. If the sample is analysed immediately on-site, then it is not necessary to dechlorinate the water.

Dechlorination

The dechlorination tablets are pre-placed inside each sample collection bag.

Simply pull off the sterile seal, pour the sample from the sampling cup into the sample collection bag, seal and shake thoroughly to ensure thorough dissolving of the tablet.

The tablet will have no adverse effect on unchlorinated water.

Tap Samples

Turn on tap and allow water to flow for about 2 minutes, to flush the sampling pipeline.

Take an initial sample and carry out any appropriate on-site physical & chemical tests e.g. Chlorine residual Turbidity, temperature & pH.

Take any other required physical and chemical samples.

Then disinfect the tap, which can be carried out in the two following ways:

Chemical Disinfection of Tap

Turn off the tap.

Squirt the inside of the tap with concentrated sodium hypochlorite solution (e.g. bleach), with a wash bottle.

Leave for 3 minutes to disinfect fully.

Flush the tap until all the bleach has been washed off – this can be verified by taking further residual chlorine tests.

Fill the prepared bacteriological sampling bottle (leave a slight air gap at the top) and seal the lid tightly.

Heat Disinfection of Tap

This method is applicable for metal taps, but not for any plastic taps or taps with non-removable anti-splash devices.

Turn off the tap fully, and flame the closed tap with a small Propane or Butane burner; cease flaming if/when any steam issues from the tap.

Flush the tap until the water cools to its original temperature.

Fill the prepared bacteriological sampling bottle (leave a slight air gap at the top) and seal the lid tightly.

Dip Sample

Sterilise the sampling cup by igniting 1ml of methanol/alcohol in the cup.

Allow to cool.

Rinse the alcohol from the cup with water from the sample source.

Immerse the cup into the water source to obtain the sample.

Pour into the prepared bacteriological sampling bottle (leave a slight air gap at the top) and seal the lid tightly.

SECTION 5: USE OF BACTERIOLOGICAL MEDIA

If stored correctly, the dissolved media should remain stable for 6-8 weeks. However, if there are any signs of contamination e.g. yellowing, cloudiness etc., discard.

Ideally, to reduce the possibility of contamination, use one bottle of media only for a 24 hour period, and use a fresh bottle on each subsequent day. However, if this is not possible, then the bottle must be resealed immediately after use. The media may be re-sterilised if required by boiling in a water bath for 15 minutes.

Clean empty media bottles thoroughly before re-use. Any residues should be washed out with hot water; cleaned with a little detergent (a small brush can be used if required); rinsed several times in clean water, dried and stored in a clean environment, with the lids lightly attached.

The MLSB solution may be applied to the pads up to 6 hours before sampling, if the pads are subsequently stored in a cool environment. This procedure can reduce the potential for contamination with excessive operations in the field.

If the MLSB powder is stored in the original sachets in dry, cool conditions it should have a shelf life of up to 5 years from the date of manufacture.

SECTION 6: ASEPTIC PROCEDURES

Aseptic procedures are of paramount importance during microbiological analysis, and extra care must be taken when outside the central laboratory, i.e. in the field.

Everything must be kept clean and sterile, particularly on the following surfaces:

Inner surface of the sampling cup

Inner surface of the graduated filter funnel

Filter membrane and absorbent pads

Upper surface of the sintered membrane support disc

Inside of the petri dishes

Support pad dispenser arm, and forceps

Before every use, rinse the filtration unit and sampling cup in clean water, and dry by using clean tissue paper.

Pour approx. 1 mL of methanol into the sampling cup and swirl. (The methanol can be stored in the plastic bottles provided).

Place the sample cup in a normal upright position, away from any flammable substances.

Using a suitable means (cigarette lighter), ignite the methanol – TAKE CARE AROUND THE NAKED FLAME.

Allow the methanol to burn for a few seconds.

Whilst the methanol is still burning, invert the filtration unit and carefully place inside the sample cup.

Wait for at least 20 minutes to ensure that the sample cup and filtration unit are sterile. Methanol burns anaerobically to form formaldehyde gas, which reaches all areas of the filtration apparatus and ensures a thorough sterilisation.

Pour any residual methanol solution away.

The above sterilisation procedures should be carried out immediately before sampling, and after the filtration of each sample.

Reusable aluminium petri dishes must be sterilized, either by immersing in boiling water, or flaming with methanol prior to use. After sterilisation, ensure that the dishes are allowed to dry thoroughly before use. Other methods of sterilisation can be employed, including autoclaving, or placing the aluminium dishes in a conventional oven at 300°C for 30 minutes. Once sterile, the petri dishes should be handled carefully to prevent subsequent recontamination.

Pads are supplied sterile, in cartridges of 100. A sterile pad dispenser is supplied for depositing the pads into the petri-dishes. It is preferable to dispense pads at the central laboratory, prior to going to the sampling point; in this way, the dispenser may be kept attached to a pad cartridge and remain clean and sterile. If it is necessary to dispense pads in the field, great care must be taken not to contaminate either the pad dispenser or the cartridge. As soon as a cartridge is finished, a new one should be attached to the dispenser. Do not leave the dispenser unattached for any length of time. If no pad dispenser is available use sterilised forceps.

7.13.1 Before handling a membrane filter with the forceps, the forceps should be flame-sterilised thusly: hold the forceps tips in a flame for at least 5 seconds, and allow to cool before handling the membrane.

SECTION 7: PROCESSING SAMPLES FOR COLIFORM ANALYSIS

All samples must be incubated within 6 hours of sampling.

Dispense a sterile absorbent pad into a sterile petri dish, and saturate the pad with prepared broth

Loosen the graduated filter funnel, and remove from the base support.

Sterilise the forceps and allow to cool. Using these forceps, place a sterile membrane onto the glass membrane support, grid side up. If the membrane tears or becomes contaminated, discard it and use a new one.



**Place the sterile membrane onto the glass membrane
using the sterilised forceps**

Lock the membrane in place by screwing the filter funnel down into position.



Pour the water sample into the filter funnel up to the 100 ml graduation.



Pour water sample up to 100 ml graduation

Connect the hand vacuum pump to the filtration unit base (ensuring the arrow is orientated to the direction of desired air flow) and pump in a controlled fashion to suck the water sample through the membrane.




Apply hand pump to pass the water through membrane

When all the water has been filtered, release the vacuum pump and use the sterile forceps to take the membrane from the filtration unit.



Use the sterile forceps to take the membrane from the filtration unit

Place the membrane on top of the pad, which has been previously saturated with the MLSB media.



Place membrane on top of MLSB saturated pad

Replace the petri-dish lid and mark the petri dish. A suitable system should be adopted to record the petri dish mark, associated with the sample number, place, date, time, etc.

Place the petri-dish into the rack. Repeat the process for all the samples and then place the filled rack into the incubator.

It is important to note that when the last sample has been processed, a resuscitation period of at least one hour (but not more than four hours) must be observed before incubating. This allows any physiologically stressed coliforms to recover before culturing.

To incubate faecal (thermotolerant) coliforms, temperature of 44.5°C should be used. For total coliform analysis, a temperature of 35°C is appropriate.


SECTION 8: COUNTING COLIFORMS AND RECORDING THE RESULT

Note the temperature that the incubator has been set for.

Following incubation, switch off the power and remove the petri dishes from the incubator.

Place the petri dishes on a flat, level surface.

Remove the lids and count all the yellow colonies, irrespective of size. Use the hand lens, if necessary. Count the colonies within a few minutes of removing from the incubator, as the colours are liable to change on cooling and standing. Ignore those colonies that are not yellow e.g. pink & transparent colonies.



Once the number of yellow colonies has been determined, and assuming that 100 mL of sample was filtered, this number of colonies equals the number of coliforms per 100 mL. Where samples were incubated at 35°C, the count is of Total Coliforms, whilst for those incubated at 44.5°C, the count is for Faecal (Thermotolerant) Coliforms.

Record the results using the record sheets provided.

SECTION 9: SELECTING THE OPTIMUM VOLUMES FOR MEMBRANE FILTRATION

The optimum volume of sample is that which will allow the most accurate quantification of bacterial colonies. This is achieved when the number of faecal (thermotolerant) coliform colonies on the membrane following incubation is between 20 and 200 colonies. If there are fewer than 10 colonies, then there exists the possibility of statistical error. Numbers greater than 200 colonies are difficult to count with the naked eye.

Potable Waters

The number of faecal coliform bacterial colonies in treated water samples should ideally be zero. Thus, the preferred sample volume is 100 mL, and a count of zero faecal coliform bacteria per 100 mL is indicative of a microbiologically safe water supply. If the count exceeds 1 faecal coliform per 100 mL, contamination is indicated. If the count exceeds 10 faecal coliforms per 100 mL, action is urgently required.


Raw Waters

For source waters and partially treated waters, including those which are ground water derived, it can be useful to adjust the sample volume in order to obtain faecal coliform counts in the optimum range 10-200. It may also be useful to process more than one quantity on the first occasion a particular water source is sampled. In such cases it is not necessary to re-sterilise the filtration equipment between different quantities of the same sample, provided that the smaller volume is processed first. Typical volumes which may be appropriate for various water types are shown in the following table. They are only guidelines; there is no substitute for experience of a given source.

SOURCE OF SAMPLE	APPROPRIATE COUNT		
	10	50	100
Lakes, Reservoirs, & Rivers &	***	**	***
Wells, boreholes, other protected	*	**	*
Water treatment plant partially	**	**	*
Water treatment plant fully treated	***		
Distribution system	***		

*** Normal Volume or First Choice ** Likely Volume

* Possible Volume



Section 5

AquaSafe

Turbidity Tube

Instruction Manual

SECTION 1: ASSEMBLY OF TURBIDITY TUBE



Carefully remove the two halves of the turbidity tube from their position, in the foam recess at the front of the AquaSafe Kit.

Align the two halves of the turbidity tube so that the graduations are easily visible, then push together the grey bolting mechanism and screw together using the grey waterproof connector.



SECTION 2: ANALYSIS

Good illumination is essential for accurate use of the turbidity tube.



Look through the open end of the tube, at the black square on the base of the tube. This is the Trace2o Secchi marker.

Hold the tube vertically, and slowly pour the water sample to be analysed into the tube, until the moment that the Secchi marker is no longer visible from the top of the tube.

Alternatively, fill the tube, then slowly pour small portions away, until the moment that the Secchi marker becomes visible from the top of the tube.


Hold the tube vertically, and identify the water level.

The turbidity value is that marked next to the graduation line nearest the water level.

The graduations follow a logarithmic scale, with the most critical values marked on the tube, and therefore the turbidity tube can only ever be an estimate for the turbidity of the water sample.

Warning – bubbles may cause false readings.

The turbidity tube is calibrated to a person with normal (6/6) visual acuity.



Section 6

AquaSafe

Comparator

Instruction Manual for WSL Range

SECTION 1: ASSEMBLY OF COMPARATOR



Insert the comparator disc into the recess on the right side of the comparator.

SECTION 2: pH ANALYSIS:



Fill **two** cells with **10mL** of sample.

Place one in the **left** side of the comparator. This is the **blank** cell.

To the other cell add one **Phenol Red Rapid Tablet**.



Crush, cap the cell and **invert to mix**. This is the **test** cell.



Place the **test** cell in the **right side** of the comparator.



Rotate the disc until colour match is obtained. Record disc reading.

pH value = disc reading

Notes:

pH values below 6.5 always produce a yellow colouration

pH values above 8.4 always produce a red colouration

Water samples with low total alkalinity may give incorrect results

Turbid samples should be filtered prior to analysis for best colour match

SECTION 3: CHLORINE ANALYSIS:



Fill **two** cells with **10mL** of sample.

Place one in the **left** side of the comparator. This is the **blank** cell.


To the other cell add one **DPD No. 1 Tablet**.



Crush, cap the cell and **invert to mix**. This is the **test** cell.



Place the **test** cell in the **right side** of the comparator.




Rotate the disc until colour match is obtained. Record disc reading of **Free Chlorine**.

Having completed test for free chlorine -

Remove cell from **right side** of comparator, remove lid and add one **DPD No. 3 Tablet**.



Crush, cap the cell and **invert to mix.**, place the **test** cell in the **right side** of the comparator.



Rotate the disc until colour match is obtained. Record disc reading of **Total Chlorine**.

Combined Chlorine mg/L (ppm) = Total Chlorine mg/L – Free Chlorine mg/L

Notes:

Take care not to shake or aerate the sample

Store the reagents in a cool, dry place

Turbid samples should be filtered prior to analysis for best colour match

APPLICATION NOTE

T2O-AN-C60

AMMONIA (C60) METHOD

The following application note explains the procedure for the detection of Ammonia (C60) using the Comparator.

Equipment:

- Comparator
- 2 x 10ml cells
- Stirring rod
- Ammonia Nos. 1 & 2 tablets

Safety:

- Consult the safety data sheet for all of the reagents before use. Even if you have used comparator reagents before, the formulation may have changed.

Getting started:**Fig 1.0 Comparator Disc**

(Disc shown is for illustration purposes only)

**Fig 1.1 Reference cell**

Insert the comparator disc for Ammonia into the recess on the right hand side of the comparator with the numbers facing forwards as shown in fig 1.0

Fill the two vials to the 10ml mark with sample water

Place one in the left side of the comparator as shown in fig 1.1. This is the reference cell

Sample preparation:


- Add **one Ammonia No.1 tablet** straight from the foil to the water sample in the second cell.
- Crush the tablet using a clean stirring rod, until no large pieces are visible.
- Add **one Ammonia No.2 tablet** straight from the foil to the same sample.
- Crush the tablet using a clean stirring rod, until no large pieces are visible.
- Fit rubber stopper to the cell and invert several times until the tablet is dissolved.
- Leave the sample to stand for 10 minutes.

Analysis:



Fig 1.2 Test cell

Place the **test** cell in the **right side** of the comparator as shown in fig 1.2



Rotate the disc until colour match is obtained. Record disc reading that is displayed in the small viewing window to the bottom right.

Ammonia concentration is the reading on the disc.

Notes

- For best results, rinse vials thoroughly between with deionised water. Ensure that the outside of the vials are clean.
- The tablets must be added in the correct sequence. Ammonia No. 1 tablet will only dissolve fully after Ammonia No. 2 tablet is added.
- The tablets should be added direct from the foil avoiding contact with hands or surfaces.
- The temperature is critical for full sample development. Below 20°C, increase reaction time to 15 minutes.
- The tablets are unsuitable for determination of ammonia in sea water or brackish water due to precipitation of salts.
- Turbid samples should be filtered prior to analysis.

APPLICATION NOTE

T2O-AN-C260

NITRATE (C260) METHOD

The following application note explains the procedure for the detection of Nitrate (C260) using the Comparator.

Equipment:

- Comparator
- 2 x 10ml cells
- Stirring rod
- Nitrate test tube
- Nitrate test powder
- Nitrate test tablet
- Nitrite LR tablet

Safety:

- Consult the safety data sheet for all of the reagents before use. Even if you have used comparator reagents before, the formulation may have changed.

Getting started:

**Fig 1.0 Comparator Disc**

(Disc shown is for illustration purposes only)

**Fig 1.1 Reference cell**

Insert the comparator disc for Nitrate into the recess on the right hand side of the comparator with the numbers facing forwards as shown in fig 1.0

Fill the two vials to the 10ml mark with sample water

Place one in the left side of the comparator as shown in fig 1.1. This is the reference cell

Sample preparation:

- Fill the nitrate test tube with 20ml of the water sample
- Add **1 level spoon of Nitrate Test Powder** (using the spoon attached to the lid of the powder pot)
- Fit the cap and swirl vigorously for one minute
- Add **one Nitrate Test tablet** straight from the foil to the nitrate test tube
- Refit the cap and swirl vigorously for 1 minute
- Stand the tube upright, allow contents to settle, then invert gently 3 to 4 times
- Allow to stand for a further 2 minutes, then open and wipe carefully around the rim to remove solid particles
- Tip 10ml of the sample from the Nitrate test tube into the second 10ml sample cell
- Add **one Nitrite LR tablet** straight from the foil to the water sample in the second cell.
- Crush the tablet using a clean stirring rod, until no large pieces are visible.
- Fit rubber stopper to the cell and invert several times until the tablet is dissolved.
- Leave the sample to stand for 10 minutes.

Analysis:



Fig 1.2 Test cell

Place the **test** cell in the **right side** of the comparator as shown in fig 1.2

Rotate the disc until colour match is obtained.

Record disc reading that is displayed in the small viewing window to the bottom right.

Nitrate concentration is the reading on the disc.



Notes

- For best results, rinse vials thoroughly between with deionised water. Ensure that the outside of the vials are clean.
- The tablets should be added direct from the foil avoiding contact with hands or surfaces.
- Nitrite present in the sample will also react leading to a higher result. To correct, carry out a nitrite test (C270) on the sample and subtract the result.
- Nitrite concentrations above 1mg/l can be diluted up to 100x, with the result multiplied up accordingly to compensate.
- If the sample flocculates it is likely the sample contains high concentrations of Nitrate. Dilute the sample at least 10x and retest.
- Interferences may occur from the presence of the following ions : Antimony (III), Iron (III), Lead, Mercury (I), Silver, Chloroplatinate, Metavanadate, Bismuth.
- Copper(II) ions may artificially lower the result.

APPLICATION NOTE

T2O-AN-C270

NITRITE (C270) METHOD

The following application note explains the procedure for the detection of NITRITE (C270) using the Comparator.

Equipment:

- Comparator
- 2 x 10ml cells
- Stirring rod
- Nitrite LR tablet

Safety:

- Consult the safety data sheet for all of the reagents before use. Even if you have used comparator reagents before, the formulation may have changed.

Getting started:

**Fig 1.0 Comparator Disc**

(Disc shown is for illustration purposes only)

**Fig 1.1 Reference cell**

Insert the comparator disc for Nitrite into the recess on the right hand side of the comparator with the numbers facing forwards as shown in fig 1.0

Fill the two vials to the 10ml mark with sample water

Place one in the left side of the comparator as shown in fig 1.1. This is the reference cell

Sample preparation:

- Add **one Nitrite LR tablet** straight from the foil to the water sample in the second cell.
- Crush the tablet using a clean stirring rod, until no large pieces are visible.
- Fit rubber stopper to the cell and invert several times until the tablet is dissolved.
- Leave the sample to stand for 10 minutes.

Analysis:



Fig 1.2 Test cell

Place the **test** cell in the **right side** of the comparator as shown in fig 1.2



Rotate the disc until colour match is obtained. Record disc reading that is displayed in the small viewing window to the bottom right.

Nitrite concentration is the reading on the disc.

Notes

- For best results, rinse vials thoroughly between with deionised water. Ensure that the outside of the vials are clean.
- The tablets should be added direct from the foil avoiding contact with hands or surfaces.
- Interferences may occur from the presence of the following ions : Antimony (III), Iron (III), Lead, Mercury (I), Silver, Chloroplatinate, Metavanadate, Bismuth.
- Copper(II) ions may artificially lower the result.

Section 7

HydroLite®

HL101

Pocket pH Tester

Instruction Manual

Introduction:


The HydroLite® HL101 is a pocket-sized pH and temperature tester. This manual provides instructions on use, with a step-by-step operating guide, as well as care and maintenance instructions.

Components:




- HydroLite® HL101 pocket pH/temperature tester
- pH Buffer Solutions (pH4.01/7.00)

Keypad:

The HL101 pocket pH/temperature tester uses a simple three-button membrane keypad, with graphical symbols to describe the function of each key.




Map of keys and functions:

KEY	FUNCTION
On/Off Hold 	<ul style="list-style-type: none"> • Power the unit ON/OFF • Freezes the currently displayed value for recording; press the key again to resume measuring. • When in calibration mode, exits calibration and returns to measurement mode.
Cal 	<ul style="list-style-type: none"> • Press the key to enter the calibration mode. • Press and hold the key to enter the setup menu. • In the setup mode, press the key to select default options.
Enter 	<ul style="list-style-type: none"> • Confirms the calibration or selected option.

Display:

The unit is equipped with a clear LCD display that is used to show measured value, units, mode indicators and function indicators.




Map of mode indicators and descriptions:


Mode indicator	Description
MEAS	Indicates unit is in the measurement mode
CAL	Indicates unit is in the calibration mode
SETUP	Indicates unit is in SETUP mode
HOLD	Indicates the measured value has been frozen
ATC	Indicates temperature compensation is enabled

Prior to Use:

- Remove the electrode protective cap from the unit.



- If the membrane on the electrode dries out, soak the electrode in 3M KCL solution or tap water for at least 15 minutes. DO NOT use distilled or deionised water, as this will shorten the life of sensor.



Power On/Off:

- Press ON/OFF key to turn on the unit, the display shows measured value.
- Press and hold the ON/OFF key for 5 seconds, the unit will turn off.
- If no key is pressed for 8 minutes, the unit will automatically turn off to conserve power.

Note: The auto-off function can be disabled if required – see below for further details

Setup Menu:

The HL101 pocket pH/temperature tester includes a comprehensive setup menu with customisable options to suit user measurement requirements.


Menu item	Description	Available options	Details	DEFAULT?
bUF	pH Buffer	USA	USA Standard (pH4.01/7.00/10.01)	•
		NIST	NIST Standard (pH4.01/6.86/9.18)	
CAL	Calibration Points	1	1 point	
		2	2 points	•
		3	3 points	
UNIT	Temperature Unit	°C	Degrees Celsius	•
		°F	Degrees Fahrenheit	
°C	Temperature Calibration	CAL	Enters the temperature calibration mode	
HOLD	Auto-Hold	YES	Automatically freezes a stable reading	
		NO	Disable	•
OFF	Auto-Off	YES	Automatically turn off the unit	•
		NO	Disable	
rst	Reset	YES	Restore factory settings	
		NO	Disable	•

Changing the default parameters:


1. Press and hold the CAL key for 3 seconds to enter the setup menu; the unit goes to buffer standard selection mode, the display shows "USA/BUF" (USA standard).




2. Press CAL key to select the USA or NIST standard for pH buffers. Press ENTER key to confirm; the unit goes into calibration point selection mode, the display shows "2/CAL" (2 points calibration).



3. Press CAL key to select the number of calibration points (1, 2 or 3 points). Press ENTER key to confirm; the unit goes into temperature unit selection mode, the display shows "°C/UNIT".



4. Press CAL key to select the temperature units (°C or °F). Press ENTER key to confirm, the unit goes into temperature calibration mode, the display shows "°C/CAL".



5. If you do not want to perform temperature calibration, press ENTER key; the unit goes to next option.

6. To perform temperature calibration, press CAL key to enter the temperature calibration mode, the display shows current temperature reading.



7. Press CAL key to set temperature value. Press ENTER key to confirm; the display shows "YES/HOLD" indicating that the auto-hold function is enabled.



If the auto-hold function is enabled, the unit will automatically sense a stable endpoint reading and freeze the reading in the display. If you disable this function, the unit allows user to freeze the reading in the display manually by pressing the HOLD key.

8. Press CAL key to enable or disable the auto-hold function. Press ENTER key to confirm; the display shows "YES/OFF" indicating that the auto-off function is enabled.



When the auto-off function is enabled, if no key is pressed for 8 minutes, the unit will automatically turn off to conserve power.

9. Press CAL key to enable or disable the auto-off function. Press ENTER key to confirm; the display shows "NO/RST" indicating the current status of the reset function.



WARNING:

The Reset function will restore the unit back to factory default settings; all calibration values and selected parameters will be reset.

10. Press CAL key to enable or disable the reset function. Press ENTER key to confirm; the unit returns to measurement mode.

EXIT THE SETUP MENU:

During the setup mode, to exit the setup menu, press ON/OFF key; the unit will return to measurement mode immediately.

pH Calibration:

The HL101 pocket pH/temperature tester allows up to 3 point calibration. We recommend that you perform at least a 2 point calibration for best accuracy. The unit automatically recognises and calibrates to the following standard buffer values.


- USA Standard Buffer Options: pH 4.01, 7.00, 10.01
- NIST Standard Buffer Options: pH 4.01, 6.86, 9.18

Single point calibration should only be carried out with pH 7.00 or pH 6.86, otherwise the calibration will not be accepted by the unit.


The unit must be calibrated prior to first use or whenever the electrode is replaced. To ensure optimum accuracy, regular calibration is recommended. Do not reuse calibration solution after calibration, contaminants in solution will affect the calibration and eventually the accuracy of the measurement.

SINGLE POINT CALIBRATION:

- 1.1 Ensure that 1 point calibration is selected in the setup menu.
- 1.2 Rinse the pH electrode with distilled water. Press CAL key; the unit shows "pH7.00/CAL1" or "pH6.86/CAL1".




- 1.3 Immerse the pH electrode in the pH7.00 buffer solution; the end of the sensor must be completely submerged in the calibration solution. Stir the solution gently.
- 1.4 Press ENTER key to confirm. Wait for the measured value to stabilise; the display shows "END". Single point calibration is completed.




2 POINT CALIBRATION:

- 2.1 Ensure that 2 point calibration is selected in the setup menu.
- 2.2 Repeat steps 1.2 to 1.4 above. When the first calibration point is completed, the display will show "CAL2". The unit prompts you to continue with second point calibration.




- 2.3 Rinse the pH electrode with distilled water. Immerse the electrode in the pH 4.01 buffer solution. Stir the solution gently.
- 2.4 Press ENTER key to confirm. Wait for the measured value to stabilize; the display shows electrode slope and "END". 2 point calibration is completed.




3 POINT CALIBRATION (not required for normal operation):


- 3.1 Ensure that 3 point calibration is selected in the setup menu.
- 3.2 Repeat steps 1.2 to 1.4 above. When the first calibration point is completed, the display will show "pH4.01/CAL2". The unit prompts you to continue with second point calibration.



- 3.3 Rinse the pH electrode with distilled water. Immerse the electrode in the pH 4.01 buffer solution. Stir the solution gently.
- 3.4 Press ENTER key to confirm. Wait for the measured value to stabilise; the display shows electrode slope and "pH10.01/CAL3".



- 3.5 Rinse the pH electrode with distilled water again. Immerse the electrode into the pH10.01 (or pH9.18) buffer solution. Stir the solution gently.
- 3.6 Press ENTER key to confirm. Wait for the measured value to stabilise; the display shows electrode slope and "END". Calibration is completed.



EXIT THE CALIBRATION:

During the calibration process, if you want to exit calibration, press ON/OFF key, the unit will return to measurement mode immediately.

Temperature Calibration:

During the measurement, if the temperature reading displayed differs from that of an accurate thermometer, you need to calibrate the unit.

1. Press and hold the CAL key for 3 seconds to enter setup menu.
2. Press ENTER key until unit shows the "°C/CAL" or "°F/CAL".



3. Press CAL key to enter the temperature calibration mode.
4. Press CAL key again to set temperature value (Resolution: 0.5°C).




5. Press ENTER key to confirm, the display shows next option.
6. Press ON/OFF key, the unit returns to measurement mode. Calibration is completed.

pH Measurement:

Rinse the pH electrode thoroughly with distilled water. Immerse the electrode into the sample solution, stir the solution gently. Wait for the reading to stabilise; record the measured value as displayed.

Hold Function:

The HL101 pocket pH/temperature tester contains two data hold modes. When the Auto-Hold function is enabled, the unit will automatically sense a stable endpoint reading and freeze it; the "HOLD" indicator appears on the display. If the Auto-Hold function is disabled, press HOLD key, the unit will immediately freeze the currently displayed value. Press the HOLD key again to resume measuring.




Electrode Care and Maintenance:

- For best results, always keep the pH electrode membrane wet.
- Ensure that the electrode is thoroughly washed with distilled water after each use.
- During extended periods of disuse, store the sensor with electrode storage solution. DO NOT use deionised or distilled water.

Electrode Replacement:

When the unit fails to calibrate or gives fluctuating readings for calibration standards, you need to replace the electrode module.

1. Twist the electrode collar counter clockwise, pull the old electrode module away from the unit.



2. Align the slot on the new electrode module, gently push the module into the unit.



3. Twist the electrode collar clockwise until it is tight. Installation is completed.

Replacing the Batteries:

If the battery indicator disappears during the use, the batteries require replacing.

1. Twist the electrode collar counter clockwise, pull the electrode module out from the unit.
2. Insert two "AAA" batteries into the battery compartment (note polarity).



3. Align the slot on the electrode module, push the electrode into the unit.
4. Twist the electrode collar clockwise until it is tight.

Troubleshooting

LCD DISPLAY	CAUSE	CORRECTIVE ACTION
---	Electrode dried out	Soak the electrode in 3M KCL solution or tap water for 10 minutes
	Measured value is out of range	Check whether the electrode membrane is clogged, dirty or broken
Err	Incorrect pH buffer solutions	Use fresh pH buffer solutions for calibration
	Electrode is broken	Replace the pH electrode module

Specifications

pH	Model	HL101
	Range	-1.00~15.00pH
	Accuracy	±0.01pH
	Resolution	0.01pH
	Calibration Points	1 to 3 points, USA (pH4.01/7.00/10.01) or NIST (pH4.01/6.86/9.18)
	Temperature Compensation	0~60°C, 32~140°F, Automatic
Temperature	Range	0~60°C, 32~140°F
	Accuracy	±1°C
	Resolution	0.1°C
	Calibration Range	Measured value ±10°C
Others	Hold Function	Manual or Automatic
	Power Off	Manual or Automatic (8 minutes after last key pressed)
	Sensor Type	Standard pH Electrode (Order Code: HL101ELEC)
	Power Requirements	2×1.5V "AAA" Batteries
	Dimensions	185(L)×40(Dia.)mm
	Weight	100g



Trace2o® Limited
The Technology Centre
Wagtech Court
Station Road
Thatcham, Berkshire, RG19 4HZ
Tel.: +44 (0)1635 566772
Fax: +44 (0)1635 873509
sales@Trace2o.com
www.Trace2o.com

The information in this document is subject to change without notice. Copyright © Trace2o® Limited, 2016. All rights reserved. Trace2o® and HydroLite® are registered trademarks of Trace2o® Ltd.

Section 8

HydroLite®

HL102

Pocket EC Tester

Instruction Manual

Introduction:


The HydroLite® HL102 is a pocket-sized conductivity (EC), TDS and temperature tester. This manual provides instructions on use, with a step-by-step operating guide, as well as care and maintenance instructions.

Components:




- HydroLite® HL102 pocket EC/TDS/temperature tester
- EC Calibration Solution (1413µS)

Keypad:

The HL102 pocket EC/TDS/temperature tester uses a simple three-button membrane keypad, with graphical symbols to describe the function of each key.




Map of keys and functions:

KEY	FUNCTION
On/Off Hold 	<ul style="list-style-type: none"> • Powers the unit ON/OFF • Freezes the currently displayed value for recording; press the key again to resume measuring. • When in calibration mode, exits calibration and returns to measurement mode.
Cal 	<ul style="list-style-type: none"> • Press the key to enter the calibration mode. • Press and hold the key to enter setup menu. • In the calibration mode, press the key to set calibration values. • In the setup mode, press the key to select default option.
Enter 	<ul style="list-style-type: none"> • Confirms the calibration or selected option. • Toggles between conductivity, TDS and salinity measurement modes.

Display:

The unit is equipped with a clear LCD display that is used to show measured value, mode indicators and function indicators.




Map of mode indicators and descriptions:


Mode indicator	Description
MEAS	Indicates unit is in the measurement mode
CAL	Indicates unit is in the calibration mode
SETUP	Indicates unit is in SETUP mode
HOLD	Indicates the measured value has been frozen
ATC	Indicates temperature compensation is enabled

Prior to Use:

- Remove the electrode protective cap from unit.



- Soak the electrode for a few minutes in tap water to remove dirt and oil stains on the electrode.






Power On/Off:

- Press ON/OFF key to turn on the unit, the display shows measured value.
- Press and hold the ON/OFF key for 5 seconds, the unit will turn off.
- If no key is pressed for 8 minutes, the unit will automatically turn off to conserve power.

Note: The auto-off function can be disabled if required – see below for further details

Setup Menu:

The HL102 pocket EC/TDS/temperature tester includes a comprehensive setup menu with customisable options to suit user measurement requirements.

PARAMETER	DESCRIPTION	OPTIONS	DESCRIPTION	DEFAULT
CAL	Calibration Point		1 point	•
			2 points	
			3 points	
TDS	TDS Factor	0.5	Setting Range: 0.1 to 1.0	0.5
UNIT	Temperature Unit	°C	Degrees Celsius	•
		°F	Degrees Fahrenheit	
°C	Temperature Calibration	CAL	Enters the temperature calibration mode	
HOLD	Auto-Hold	YES	Automatically freezes a stable reading	
		NO	Disable	•
OFF	Auto-Off	YES	Automatically turn off the unit	•
		NO	Disable	
RSt	Reset	YES	Restore factory settings	
		NO	Disable	•

SETTING THE DEFAULT PARAMETERS:

1. Press and hold the CAL key for 3 seconds to enter setup menu; the display shows currently selected number of calibration points.




2. Press CAL key to select 1 or 2 or 3 point calibration. Press ENTER key to confirm; the unit goes to TDS factor setting mode, the display shows "0.5/TDS".



3. Press CAL key to set the TDS conversion factor. Press ENTER key to confirm; the unit goes to temperature unit selection mode, the display shows "°C/UNIT".



4. Press CAL key to select the temperature unit (°C or °F). Press ENTER key to confirm; the unit goes to temperature calibration mode, the display shows "°C/CAL".



5. If you do not want to perform temperature calibration, press ENTER key; the unit goes to next option.
6. To perform temperature calibration, press CAL key to enter the temperature calibration mode, the display shows current temperature reading.



7. Press CAL key to set temperature value. Press ENTER key to confirm; the display shows "YES/HOLD" indicating that the auto-hold function is enabled.



If the auto-hold function is enabled, the unit will automatically sense a stable endpoint reading and freeze the reading in the display. If you disable this function, the unit allows user to freeze the reading in the display manually by pressing the HOLD key.

8. Press CAL key to enable or disable the auto-hold function. Press ENTER key to confirm; the display shows "YES/OFF" indicating that the auto-off function is enabled.



When the auto-off function is enabled, if no key is pressed for 8 minutes, the unit will automatically turn off to conserve power.

9. Press CAL key to enable or disable the auto-off function. Press ENTER key to confirm; the display shows "NO/RST" indicating the current status of the reset function.



WARNING:

The Reset function will restore the unit back to factory default settings, all calibration values and selected parameters will be reset.

10. Press CAL key to enable or disable the reset function. Press ENTER key to confirm, the unit returns to measurement mode. Setting is completed.

EXIT THE SETUP MENU:

During the setup mode, to exit the setup menu, press ON/OFF key; the unit will return to measurement mode immediately.

Conductivity Calibration


The HL102 pocket EC/TDS/temperature tester allows up to 3 point calibration in the conductivity mode. To ensure higher accuracy, we recommend that you perform a 3 point calibration, or select a calibration standard concentration close to the sample value you are measuring. For typical environmental freshwater samples, 1413 μ S/cm should be sufficient. The unit will automatically detect these conductivity standard solutions and prompt the user to calibrate the meter. When the calibration is done, all new calibration values will automatically override existing data. The following table shows acceptable conductivity ranges of calibration solution for each measuring range.

MEASURING RANGE	CALIBRATION SOLUTION RANGE	DEFAULT
0~20 μ S/cm	7~17 μ S/cm	10 μ S/cm
20~200 μ S/cm	70~170 μ S/cm	84 μ S/cm
200~2000 μ S/cm	700~1700 μ S/cm	1413 μ S/cm
2~20mS/cm	7~17mS/cm	12.88mS/cm

Ensure that you use fresh conductivity standard solution during the calibration. Do not reuse calibration solutions as it may be contaminated and affect the calibration and accuracy of measurement.

SINGLE POINT CALIBRATION:

- 1.1 Rinse the conductivity electrode with distilled water, then rinse with a small amount of calibration solution.
- 1.2 Press CAL key; the unit enters calibration mode.




- 1.3 Immerse the conductivity electrode in the calibration solution; the unit immediately displays current calibration standard (e.g., 1413 μ S/cm).
- 1.4 Press ENTER key; the default calibration value begins flashing.



- 1.5 Press CAL key to set each digit, press ENTER key to confirm. When the setting is done, ensure the displayed value matches the chosen calibration standard.




- 1.6 Press ENTER key to start the calibration. Wait for the reading to stabilise; the display shows "END". The unit returns to measurement mode automatically. Single point calibration is completed.



2 POINT CALIBRATION:

- 2.1 Ensure that 2 point calibration is selected in the setup menu.
2.2 Repeat steps 1.2 to 1.6 above; when the first calibration point is done, the display will show "CAL2". The unit prompts you to continue with second point calibration.



- 2.3 Immerse the conductivity electrode in the calibration solution, the unit automatically shows current calibration standard (e.g., 84µS/cm).



- 2.4 Press ENTER key; the default calibration value begins flashing.
2.5 Press CAL key to set each digit, press ENTER key to confirm. When the setting is done, ensure the displayed value matches the chosen calibration standard.
2.6 Press ENTER key to start the calibration. Wait for the reading to stabilise; the display shows "END". The unit returns to measurement mode. Second point calibration is completed.

3 POINT CALIBRATION:


- 3.1 If 3 point calibration is selected in the setup menu, when second calibration point is done, the display will immediately show "CAL3". The unit prompts to continue with third point calibration.
- 3.2 Repeat the above steps until the display shows "END", the unit returns to measurement mode. Calibration is completed.

Note: Performing the conductivity calibration will simultaneously calibrate the corresponding TDS and salinity values.


Temperature Calibration

During the measurement, if the temperature reading displayed differs from that of an accurate thermometer, you need to calibrate the unit.

1. Press and hold the CAL key for 3 seconds to enter setup menu.
2. Press ENTER key until unit shows the "°C/CAL" or "°F/CAL".



3. Press CAL key to enter the temperature calibration mode.
4. Press CAL key again to set temperature value (Resolution: 0.5°C).




5. Press ENTER key to confirm, the display shows next option.
6. Press ON/OFF key, the unit returns to measurement mode. Calibration is completed.

EXIT THE CALIBRATION:


During the calibration process, if you want to exit the calibration, press ON/OFF key, the unit will return to measurement mode.

Switching Measurement Mode:


1. In the conductivity mode, press ENTER key until the display shows “TDS”. The unit enters TDS measurement mode.



2. Press ENTER key, the unit enters the salinity measurement mode.



3. Press ENTER key again, the unit returns to conductivity measurement mode.




Measurement:

Rinse the conductivity electrode with distilled water. Immerse the electrode into the sample solution. Stir the solution gently. Wait for the reading to stabilise; record the measured value as displayed.

Hold Function:

The HL102 pocket EC/TDS/temperature tester contains two data hold modes. When the Auto-Hold function is enabled, the unit will automatically sense a stable endpoint reading and freeze it; the “HOLD” indicator appears on the display. If the Auto-Hold function is disabled, press HOLD key, the unit will immediately freeze the currently displayed value. Press the HOLD key again to resume measuring.




Electrode Care and Maintenance

- After each use, the conductivity electrode should be rinsed thoroughly in deionised water.
- If there is a build-up of solids inside the measurement area of the cell, these should be removed very carefully with a cotton bud soaked in solvent, taking care not to touch the metal parts of the inner cell.

Electrode Replacement:

When the unit fails to calibrate or gives fluctuating readings for calibration standards, you need to replace the electrode module.

1. Twist the electrode collar counter clockwise, pull the old electrode module away from the unit.



2. Align the slot on the new electrode module, gently push the module into the unit.



3. Twist the electrode collar clockwise until it is tight. Installation is completed.

Replacing the Batteries:

If the battery indicator disappears during the use, the batteries require replacing.

1. Twist the electrode collar counter clockwise, pull the electrode module out from the unit.
2. Insert two “AAA” batteries into the battery compartment (note polarity).



3. Align the slot on the electrode module, push the electrode into the unit.
4. Twist the electrode collar clockwise until it is tight.

Troubleshooting

LCD DISPLAY	CAUSE	CORRECTIVE ACTION
■■■	Electrode dried out	Soak the conductivity electrode in tap water for 10 minutes
	Measured value is out of range	Check the electrode whether clogged, dirty or broken
Err	Incorrect calibration solutions	Using the fresh conductivity standard solutions for calibration
	Setting value does not match calibration solution	Reset the calibration value
	Electrode is broken	Replace the electrode module

Specifications

Conductivity	Model	HL102
	Range	0~20.00, 200.0, 2000μS/cm, 20.00mS/cm
	Accuracy	±1% F.S
	Resolution	0.01, 0.1, 1
	Calibration Points	1 to 3 points
	Calibration Solutions	10μS/cm, 84μS/cm, 1413μS/cm, 12.88mS/cm
TDS	Range	0~10ppt (Max. 20ppt, depending on factor setting)
	Accuracy	±1% F.S
	TDS Factor	0.1~1.0 (Default 0.5)
Salinity	Range	0~10ppt
	Accuracy	±1% F.S
	Resolution	0.01ppt
Temperature	Range	0~60°C, 32~140°F
	Accuracy	±1°C
	Resolution	0.1°C, 0.1°F

	Calibration Points	1 point
	Calibration Range	Measured value $\pm 10^{\circ}\text{C}$
Others	Temperature Compensation	0~60°C, 32~140°F
	Temperature Coefficient	2%/°C
	Cell Constant	K=1
	Normalization Temperature	25°C
	Hold Function	Manual or Automatic
	Power Off	Manual or Automatic (8 minutes after last key pressed)
	Sensor Type	Order Code: HL102ELEC
	Power Requirements	2×1.5V "AAA" Batteries
	Dimensions	185(L)×40(Dia.)mm
	Weight	100g

Calculation of TDS conversion factor:

To determine the TDS conversion factor use the following formula:

$$\text{Factor} = \frac{\text{Actual TDS}}{\text{Actual Conductivity @ } 25^{\circ}\text{C}}$$

Where:

Actual TDS: value from the high purity water and precisely weighed NaCl or KCL reagent.

Actual Conductivity: the measured conductivity value.

For example:

Dissolve 64 grams of potassium chloride reagent in 1L distilled water. If its conductivity value is 100mS/cm, then TDS conversion factor is 0.64.



Trace2o® Limited
The Technology Centre
Wagtech Court
Station Road
Thatcham, Berkshire, RG19 4HZ
Tel.: +44 (0)1635 566772
Fax: +44 (0)1635 873509
sales@Trace2o.com
www.Trace2o.com

The information in this document is subject to change without notice. Copyright © Trace2o® Limited, 2016. All rights reserved. Trace2o® and HydroLite® are registered trademarks of Trace2o® Ltd.